



---


**BIO-RAD**

**Professional Development**

## Novel Research Starts with... *GAPDH*



## Instructors



### **Stan Hitomi**

Coordinator – Math & Science

Principal – Alamo School

San Ramon Valley Unified School District

Danville, CA

### **Kirk Brown**

Lead Instructor, Edward Teller Education Center

Science Chair, Tracy High School

and Delta College, Tracy, CA

### **Bio-Rad Curriculum and Training Specialists: Sherri Andrews, Ph.D.**

[sherri\\_andrews@bio-rad.com](mailto:sherri_andrews@bio-rad.com)

### **Essy Levy, M.Sc.**

[essy\\_levy@bio-rad.com](mailto:essy_levy@bio-rad.com)

### **Leigh Brown, M.A.**


[leigh\\_brown@bio-rad.com](mailto:leigh_brown@bio-rad.com)

## Why Teach

### Novel Research Using *GAPDH* PCR?


- **Students conduct inquiry-based experiments**
- **Students understand research is a process rather than a single experiment**
- **Students learn laboratory skills and techniques commonly used in research**
- **Students guide the research process and make decisions about their next steps**
- **Students formulate scientific explanations using data, logic, and evidence**
- **Students learn from failures and unexpected results**

## Workshop Timeline



- **Introduction**
- **Preparation of Initial PCR Reactions**
- **Exonuclease Treatment**
- **Preparation of Nested PCR Reactions**
- **Analyze PCR results**

# Laboratory Overview




## Benefits of using plants



- **Large number of species**
- **Lots of diversity**
- **Phylogenetic approaches**
- **Avoid ethical concerns associated with animals**
- **No pre-approval**



## What is a Housekeeping Gene?



**Highly conserved genes that must be continually expressed in all tissues of organisms to maintain essential cellular functions.**


### **Examples:**

- ***GAPDH***
- ***Cytochrome C***
- ***ATPase***
- ***β-actin***




# Glyceraldehyde 3-Phosphate Dehydrogenase (GAPDH)

- Enzyme of glycolysis
- Structure and reaction mechanism well-studied
- Multitude of sequences
- Highly conserved



# Gene Families




Enzyme Function	GAPDH Protein Subunit	EC Designation*	Arabidopsis Gene
NAD <sup>+</sup> -dependent GAPDH in cytosol	GAPC	EC 1.1.1.12	GAPC
	GAPC-2		GAPC-2
NAD <sup>+</sup> -dependent GAPDH in plastids	GAPCP	EC 1.1.1.12	GAPCP
	GAPCP-2		GAPCP-2

## Extracting DNA


How do we isolate and identify the *GAPDH* gene?

- Use young, fresh plant-tissue
- DNA extraction at room temperature
- Time requirement ~30 minutes
- Does not require DNA quantification




- DNA can be extracted from the fresh plant tissue and multiple copies of the *GAPDH* gene can be amplified using PCR

## What do we need for PCR?



- **Template DNA (extracted from plants)**
- **Nucleotides (dATP, dCTP, dGTP, dTTP)**
- **Taq DNA polymerase**
- **Magnesium chloride (enzyme cofactor)**
- **Buffer, containing salt**
- **Sequence-specific primers flanking the target sequence**



## DNA sequence varies between species


## Primers need to be designed to account for species variation

		Section 85											
	(3781)	3781	3790	3800	3810	3825							
F309.10 CUT FROM F309	(1705)	ACCTTTTCATTTGTCCT-GT	-----G	CAGGGT	GGTGC	AA	GAA	GGT					
T8K14.5 cut from T8K14 or AC007202	(1629)	ATATT-CATCGTTT-TT	-----A	TAGGGC	GGTGC	CAA	GAA	AGT					
DQ075672 oxalis stricta	(112)	ATATAAATTTTGAATT	-----T	CAGGGT	GGTGC	AAA	GAA	AGT					
Lime AB 2 T7 minus pTZ rc	(143)	-TCTTGACTTTTCT-AT	-----A	CAGGGT	GGTGC	AAA	GAA	AGT					
OJ1116_A06.36 cut from AP004113-rice GAPDH	(2924)	-TTATGCTTTGTTT	-----G	TAGGGC	GGTGC	AA	GAA	AGT					
GAPDH cut from CAN272042-capsicum annuum	(3636)	CTGTTTGAATTGT	-----G	CAGGGT	GGTGC	AAA	GAA	AGT					
GAPDH cut from ZMGP1-maize	(1627)	TACTTTCTTTTGAATG	-----A	CAGGGT	GGTGC	CAA	GAA	AGT					
AB110609 lycopersicon esculentum	(421)	TAATTTACTTGGATGTTTGTCTGT	-----C	AGGGT	GGTGC	AAA	GAA	AGT					
Consensus (3781)		T TTTCTTGT T TT		CAGGGTGGTGCC	AAA	GAA	AGT						

		Section 86									
	(3826)	3826	3840	3850	3860	3870					
F309.10 CUT FROM F309	(1742)	CATCATTTCTGCACCTTCAGCGGACGGGCCCATGTTTGT	TGT	TGG							
T8K14.5 cut from T8K14 or AC007202	(1665)	TATAATTTCTGCACCTTCTGOTGACGCCACCATGTTTGT	TGT	TGG							
DQ075672 oxalis stricta	(150)	GATCATTTTCAGCTCCATCAGCTGATGCGCCATGTTTGT	G	TGG							
Lime AB 2 T7 minus pTZ rc	(179)	GGTAATATCAGCTCCATCAGCTGATGCACCAATGTTTGT	G	TAGG							
OJ1116_A06.36 cut from AP004113-rice GAPDH	(2958)	TGTGATATCTGCTCCATCAGCAGATGCTCCAATGTTTGT	CT	TGG							
GAPDH cut from CAN272042-capsicum annuum	(3669)	CGTAATCTCAGCCCCATCAGCTGATGCACCAATGTTTGT	G	TGGG							
GAPDH cut from ZMGP1-maize	(1665)	TGTTATCTCTGCACCTAGCAAAACGCCACCATGTTTGT	TGT	TGG							
AB110609 lycopersicon esculentum	(466)	TGTGATCTCTGCTCCATAGCAAAAGATGCTCCCATGTTTGT	TGT	TGGG							
Consensus (3826)		TGT AT TCTGCTCCATCAGCTGATGCACC	ATGTTTGT	TGT	TGG						

- What if you don't know the exact DNA target sequence?
- How do you design primers?



## Genetic code Degeneracy



- **Multiple codons code for one amino acid**
- **Variations at the DNA level may not be reflected at the **protein** level**

### Amino acids and the DNA codons for each.

Ala (A)	Arg (R)	Asp (D)	Asn (N)	Cys (C)	Gln (E)	Glu (Q)	Gly (G)	His (H)	Ile (I)
GCA	CGA	GAC	AAC	TGC	CAA	GAA	GGA	CAC	ATA
GCC	CGC	GAT	AAT	TGT	CAG	GAG	GGC	CAG	ATC
GCG	CGG						GGG		ATT
GCT	CGT						GGT		
	AGA								
	AGG								
Leu (L)	Lys (K)	Met (M)	Phe (F)	Pro (P)	Ser (S)	Thr (T)	Trp (W)	Tyr (Y)	Val (V)
CTA	AAA	ATG	TTC	CCA	TCA	ACA	TGG	TAC	GTA
CTC	AAG		TTT	CCC	TCC	ACC		TAT	GTC
CTG				CCG	TCG	ACG			GTG
CTT				CCT	TCT	ACT			GTT
TTA					AGC				
TTG					AGT				

**Note:** When there are multiple codons for an amino acid, the codons are very similar. For all amino acids with up to four codons, only the third base differs between codons (for example, the four codons for valine, which all begin with GT). Also three codons code for a stop signal, these are TAG, TGA, TAA.

## Variation in the Genetic Code



**Introns are less conserved since they do not code for protein**

**The GAPDH enzyme (protein) is highly conserved but there are variations at the DNA level**



**Conservative substitution-doesn't change protein properties**



## Degenerate primers are used to account for sequence variation

Plant	Gene	Accession number	Sequence
<i>Arabidopsis</i>	<i>GAPC1</i>	AT3G04120	GACTACGTTGTTGAGTCTACTGG
<i>Arabidopsis</i>	<i>GAPC2</i>	AT1G13440	GACTTTGTTGTTGAGTCTACTGG
<i>Arabidopsis</i>	<i>GAPCP1</i>	AT1G79530	GATTATGTTGTTGAGTCTTCCGG
<i>Arabidopsis</i>	<i>GAPCP2</i>	AT1G16300	GAGTATGTTGTTGAGTCTTCCAGG
Pepper	<i>GAPCP</i>	CAN272042	GATTATGTTGTTGAATCTTCTGG
Liverwort	<i>GAPC</i>	AJ246023	GAGTACGTCGTCCGAGTCTACCGG
Corn	<i>GAPC1</i>	ZMGPC1	GAGTACGTCGTGGAGTCCACCGG
Corn	<i>GAPC2</i>	U45855	GAGTATGTCGTGGAGTCCACCGG
Corn	<i>GAPC3</i>	U45856	GAATATGTTGTTGAGTCTACTGG
Corn	<i>GAPC4</i>	X73152	GAATATGTTGTTGAGTCTACTGG
Pea	<i>GAPC1</i>	L07500	GATATCATTGTTGAGTCTACTGG
Wheat	<i>GAPC</i>	EF592180	GAGTACGTTGTTGAGTCCACCGG
Rye Grass	<i>GAPC3</i>	EF463063	GACTACGTTGTTGAGTCCACTGG
Tobacco	<i>GAPC</i>	AJ133422	GATTACATTGTGGAGTCCGACTGG
Tobacco	<i>GAPDH*</i>	DQ682459	GATTTGTTGTGGGAATCCACTGG
Carrot	<i>GAPDH*</i>	AY491512	GAGTACATTGTGGAGTCCACTGG
Blue Gem	<i>GAPDH*</i>	X78307	GAGTACGTCGTTGAGTCCGACTGG
Tomato	<i>GAPDH*</i>	AB110609	GACTTCGTTGTTGAATCAACCGG
Snapdragon	<i>GAPDH*</i>	X59517	GAGTATATTGTGGAGTCCACTGG
<b>INITIAL FORWARD PRIMER</b>			GABTATGTTGTTGARTCTTCWGG
Position of base			123456789

## Primers are designed using the consensus sequence

Plant	Accession Number	Sequence
<i>GAPDH</i>		
Tobacco	DQ682459	GATTTCGTTGTGGAATCCACTGG
Carrot	AY491512	GAGTACATTGTGGAGTCCACTGG
Blue gem	X78307	GAGTACGTCGTTGAGTCGACTGG
Tomato	AB110609	GACTTCGTTGTTGAATCAACGG
Snapdragon	X59517	GAGTATATTGTGGAGTCCACTGG

Consensus sequence GABTATGTTGTTGARCTCTTCWGG

Primer set:

GA(GTC)TATGTTGTTGA(GA)TCTTC(AT)GG

Yield: 12 primers

Reverse primers are designed in the same fashion

## Degenerate primers have optional bases in specified positions

IUB Code	Bases	Derivation of IUB Code
N	A/G/C/T	any
K	G/T	keto
S	G/C	strong
Y	T/C	pyrimidine
M	A/C	amino
W	A/T	weak
R	G/A	purine
B	G/T/C	not A
D	G/A/T	not C
H	A/C/T	not G
V	G/C/A	not T

**To increase the probability that the primer will anneal to the target DNA, variable bases are designed into the primer.**

## Multiple oligos comprise the forward primer

☆  
 GAGTATGTTGTTGA (GA) TCTTC (AT) GG } **Position 3**  
 GATTATGTTGTTGA (GA) TCTTC (AT) GG } **has 3 bases**  
 GACTATGTTGTTGA (GA) TCTTC (AT) GG }

GA (GTC) TATGTTGTTGA ☆ GCCTTC (AT) GG } **Position 15**  
 GA (GTC) TATGTTGTTGA ATCTTC (AT) GG } **has 2 bases**

GA (GTC) TATGTTGTTGA (GA) TCTTC ☆ AGG } **Position 21**  
 GA (GTC) TATGTTGTTGA (GA) TCTTC TGG } **has 2 bases**

$$3 \times 2 \times 2 = 12$$

**different oligonucleotides comprising the  
forward primer**

## DNA Isolation and Amplification




To identify differences in *GAPDH* code we must isolate plant DNA and amplify the gene of interest using PCR first with **degenerate primers** (primers that account for variation in the DNA code)

A second PCR reaction (**Nested PCR**) is necessary to amplify the region which contains **one** of the *GAPDH* gene sequences (second set of primers are **nested** inside the initial PCR product sequence)

Biotechnology Explorer PCR primers are color-coded

# Setting up initial PCR Reactions



## GAPDH PCR – Quick Guide

### Stage 2A

1. Label five PCR tubes with your initials and the following labels:  
1 – Negative control (sterile water)  
2 – *Arabidopsis* gDNA  
3 – Positive control pGAP plasmid  
4 – Plant 1 gDNA (optional)  
5 – Plant 2 gDNA (optional)




2. Place each PCR tube into a tube adaptor, cap each tube. Place the adaptor tube with PCR tube on ice.




3. Flick to thoroughly mix reagent tubes and centrifuge 10 sec to force contents to bottom of tubes.



4. Pipet 20 µl of 2x blue mastermix with initial primers (2x MMIP) into each PCR tube.



5. Using a fresh tip, pipet 15 µl of sterile water to each tube.




6. Using a fresh tip every time, add 5 µl of the appropriate DNA template to each tube. Gently pip up and down to mix reagents. Recap tubes.



7. When instructed, place PCR tubes in thermal cycler.



## Setting up Initial PCR Reactions




### Protocol

1. Add 20  $\mu$ l of blue mastermix with initial primers to each PCR tube
2. Add 15  $\mu$ l of sterile water to each tube
3. Add 5  $\mu$ l of DNA template to the appropriate tube
  - Control Arabidopsis gDNA
  - Control plasmid DNA
  - Test gDNAs
  - Negative control
4. Amplify in thermal cycler (Annealing temp 52°C)

## Initial PCR Reaction

(with degenerate primers)



**PCR Animation**


[http://www.bio-rad.com/flash/07-0335/07-0335\\_PCR.html](http://www.bio-rad.com/flash/07-0335/07-0335_PCR.html)



## Results of Initial PCR Reactions

**1% agarose gel loaded with (20 µl) initial PCR samples.**

**Green bean and Lamb's ear gDNA samples generated using Nucleic Acid Extraction module.**



Lane 1- 500 bp molecular weight ruler (10 µl),

Lane 2- PCR of **control *Arabidopsis* gDNA** with initial primers (20 µl)

Lane 3- PCR of **green bean gDNA** with initial primers (20 µl)


Lane 4- PCR of **lamb's ear gDNA** with initial primers (20 µl)

Lane 5- PCR of **pGAP plasmid control** with initial primers (20 µl)


**Sometimes no amplification is observed with initial PCR with some plants, or much fainter than this gel**

## Initial PCR may result in some amplicons that are non-specific

Because of the degenerate primer used to amplify the *GAPDH* of various plant species the initial PCR may also result in some non-specific amplifications



**Nested PCR is used to amplify specifically within the amplicon (PCR product) of interest**



## Why is a Nested PCR reaction necessary?


- **Use of degenerate primers may not give you an exact match to the target sequence**

INITIAL FORWARD PRIMER	GACTATGTTGTTGAGTCTTCTGG
Arabidopsis GAPC1AT3G04120	GACTACGTTGTTGAGTCTACTGG

- **Because there are multiple primers in the mix, the primer concentration for the matching primer is lower than normal (1/12<sup>th</sup> concentration)**
- **Problems with initial PCR:**
  - inefficient
  - non-specific
- **Benefits of initial PCR:**
  - cast a wide net
  - increase the pool of specific products


## Using Nested PCR to increase your final PCR product

### Initial PCR



**DNA template:  
Genomic DNA**

### Nested PCR





**DNA template:  
Initial PCR products**

- **There is more PCR product from the nested PCR reactions since there is more specific template DNA to start from**
- **Results: intense, bold band on agarose gel**

## Exonuclease 1 treatment is needed before the second round of PCR (nested PCR) is done

- **The primers that were not incorporated into PCR product in the first reaction must be removed so that they do not amplify target DNA in the second round of PCR.**
- **Exonuclease I will be added to the PCR products**
- **The enzyme must be inactivated before proceeding to the nested PCR**

## Exo treatment




### Protocol

- 1. Add 1  $\mu$ l of exonuclease 1 enzyme to the PCR reactions**
- 2. Incubate 15 min 37°C**
- 3. Incubate 15 min 80°C**
- 4. Dilute 2  $\mu$ l Exo-treated PCR product in 98  $\mu$ l water**

**Note: Thermal cycler can be programmed for exo incubations**

**Nested PCR  
amplifies only  
regions within  
the *GAPDH*  
gene**

**Nested PCR is  
more specific**





## Setting up Nested PCR Reactions




### Protocol

- 1. Mix 20  $\mu$ l of diluted exo-treated template DNA with 20  $\mu$ l of yellow mastermix with nested primers**
- 2. For controls, mix 20  $\mu$ l of control pGAP plasmid and 20  $\mu$ l of water with 20  $\mu$ l of yellow mastermix with nested primers**
- 3. Amplify in thermal cycler (Annealing temp 46°C)**


# PCR results

**1% agarose gel loaded with 20  $\mu$ l initial PCR samples and 5  $\mu$ l nested PCR samples.**




## PCR Products

## Initial vs. Nested Reactions




<i>Arabidopsis</i> GAPC Gene	Length of PCR product (bp)	
	Initial primers	Nested primers
<i>GAPC</i>	1065	993
<i>GAPC-2</i>	1216	1145
<i>GAPCP-1</i>	1303	1231
<i>GAPCP-2</i>	1205	1133



## Cloning and Sequencing Explorer Series GAPDH PCR Module

### Results of Initial and Nested PCR




- Lane 1- 500 bp molecular weight ruler (10  $\mu$ l),
- Lane 2- *Arabidopsis* gDNA with initial primers (20  $\mu$ l)
- Lane 3- Initial *Arabidopsis* PCR product with nested primers (5  $\mu$ l)
- Lane 4- Green bean gDNA with initial primers (20  $\mu$ l)
- Lane 5- Initial Green bean PCR product with nested primers (5  $\mu$ l)
- Lane 6- Lamb's ear gDNA with initial primers (20  $\mu$ l)
- Lane 7- Initial Lamb's ear PCR product with nested primers (5  $\mu$ l)
- Lane 8- pGAP plasmid control with initial primers (20  $\mu$ l)
- Lane 9- pGAP plasmid control with nested primers (5  $\mu$ l)

**BIO-RAD**

## Cloning and Sequencing Explorer Series GAPDH PCR Module

### Results of Initial PCR



- Lane 1- 500 bp molecular weight ruler (10  $\mu$ l),
- Lane 2- Control *Arabidopsis* gDNA with initial primers (20  $\mu$ l)
- Lane 3- Green bean gDNA with initial primers (20  $\mu$ l)
- Lane 4- Lamb's ear gDNA with initial primers (20  $\mu$ l)
- Lane 5- pGAP plasmid control with initial primers (20  $\mu$ l)

**BIO-RAD**

## Webinars

- **Enzyme Kinetics — A Biofuels Case Study**
- **Real-Time PCR — What You Need To Know and Why You Should Teach It!**
- **Proteins — Where DNA Takes on Form and Function**
- **From plants to sequence: a six week college biology lab course**
- **From singleplex to multiplex: making the most out of your realtime experiments**

**[explorer.bio-rad.com](http://explorer.bio-rad.com) → Support → Webinars**