

# DISTEL

— C —

## DOSING SYSTEM

### Microbiological Efficacy Summary

	ORGANISM	TEST NORM	TEST TYPE	CONDITIONS	CONTACT TIME
VIRUSES	Bovine Viral Diarrhoea Virus (surrogate of Hepatitis C Virus)	DVV and RKI	Suspension	Dirty <sup>1</sup>	3 min.
	Modified Vaccinia Virus Ankara				
FUNGI / YEAST	<i>Candida albicans</i>	EN 13624	Suspension	Clean <sup>1</sup>	3 min.
BACTERIA	<i>Staphylococcus aureus</i>	EN 13727	Suspension	Clean <sup>1</sup>	3 min.
	<i>Pseudomonas aeruginosa</i>				
	<i>Enterococcus hirae</i>				
	Carbapenem-resistant Enterobacteriaceae (CRE) <i>Klebsiella pneumoniae</i>				
	Vancomycin-resistant Enterococci (VRE) <i>Enterococcus faecium</i>				
	<i>Escherichia coli</i>				
	Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)				
	Multi-drug resistant <i>Acinetobacter baumannii</i> (MDRAB)				
DNA / RNA	<b>TEST METHOD</b>				
	Spectrophotometric and fluorometric				

**Clean/Dirty Conditions Key:**

**Clean 1:** Testing performed in the presence of 0,3 g/L bovine albumin.

**Dirty 1:** Testing performed in the presence of 10% foetal calf serum.