



RT- qPCR > Owi Green FAST qPCR Premix

Reliability. Accuracy. Success.


We are with you to improve your research.

PCR

Polymerase (Standard, Hot Start, High Fidelity), dNTP mix

RT-qPCR

- **Owi Green FAST qPCR Premix**
- Owi Script RT UltraMix with gDNA CleanUp for qPCR
- Owi dNTP Mix 40mM (10mM dATP, dCTP, dGTP & dTTP)



ELECTROPHORESIS

Agarose, Dye DNA Stain, DNA Ladder

WESTERN BLOT


ECL, Markers

Owi Green FAST qPCR Premix

Premix 2x including Owi Green (SYBR®/FAM channel detection)


- **Universal:** compatible with all qPCR instruments ; Rox I / Rox II or without Rox
- **Performance:**
 - Sensitivity - fewer false negatives
 - Specificity - more reliable results
 - Better resolution - more accurate and sensitive results
 - Dynamic range of 6-Log
- **Rapid protocol**
- **Formats in 1 ml, 5 ml, 5 x 5 ml :** 2x Owi Green Premix, Rox I and Rox II

Excellent sensitivity



Template : 6-log dilutions of mouse gDNA. The target mouse GAPDH gene was detected by **Owi Green Fast qPCR Premix**. The results demonstrate that **Owi Green Fast qPCR Premix** can accurately amplify all 6-log dilutions between 12 – 32 Ct, showing broad range amplification ability.

Superior specificity



Template : 6-log gradient dilutions of human gDNA. The target TBP gene was detected by **Owi Green Fast qPCR Premix** and by the competitors T, V and M. The results show that the **Owi Green Fast qPCR Premix** gives more specific amplification and more accurate results.

Reference	Name	Size
OZYA008-1	Owi Green FAST qPCR Premix	1 ml / 100 reactions (20 µl)
OZYA008-5XL	Owi Green FAST qPCR Premix	5 ml / 500 reactions (20 µl)
OZYA008-25	Owi Green FAST qPCR Premix	5 x 5 ml / 2500 reactions (20 µl)

[VISIT OUR WEBSITE TO ORDER](#)

Your local vendor :