



WESTERN BLOT > Owi EColor® Prestained Protein Marker

Reliability. Accuracy. Success.


We are with you to improve your research.

PCR

Polymerase (Standard, Hot Start, High Fidelity), dNTP mix

RT-qPCR

Master Mix, dNTP mix




ELECTROPHORESIS

*Agarose, Dye DNA Stain,
DNA Ladder*

WESTERN BLOT

- **Owi EColor® Prestained Protein Marker**
- Owi ECL RevelBIOT® Pico
- Owi ECL RevelBIOT® Femto

Owi EColor® Prestained Protein Marker



6.5 kDa - 270 kDa / SDS-PAGE Tris-Glycine

- **Universal:** suitable for large and small proteins
- **Easy spotting:** tricolor (7 blue, 2 orange & 1 green bands)
- **Accurate:** high band resolution on gel
- **Quality:** premium proteins and covalent dye conjugation
- **Eco-friendly:** ecological & economical
 - 2 weeks at 25°C
 - 3 months at 4°C
 - 2 years at - 20°C

Reference	Name	Size
OZYC005-25	Owi EColor® Prestained Protein Marker	5 lanes for 10-well gel (25 µL)
OZYC005-500	Owi EColor® Prestained Protein Marker	100 lanes for 10-well gel (500 µL)

[VISIT OUR WEBSITE TO ORDER](#)