Melting Point

1. Cole-Parmer[®] Stuart[™] MP-200D **Digital Melting Point Apparatus**

Designed for maximum safety and ease of use. Select, measure, and display the temperature digitally, ensuring accuracy and avoiding the need for a glass thermometer. Test two samples simultaneously, viewing the melt via a magnifying lens with built in illumination. Set a plateau temperature for initial warm up - unattended warm up time frees up the user until the plateau temperature is reached. Heating rate continues at 2°C after the warm up, slowly reaching the melt temperature for greater accuracy.

- Easily view the melt via magnifying lens with built in illumination
- Be more productive test two samples simultaneously
- Extend the back feet for the optimal viewing angle
- 3 digit LED display
- Supplied with 100 x capillary tubes (open at one end)

Code	Pack	Price
MEL2104	Each	£754.95





PHYSICS

Т V

W